

CONTENU DU KIT





INSTI[®]
HCV Antibody Test

Introducing the World's First One Minute
HCV Antibody Test

- + Easy to use with results in one minute
- + Clear and accurate results
- + Enhances patient flow
- + Minimizes face to face contact time






HCV Antibody Test

Detects HCV genotypes 1 - 6


IgG/IgM Procedural Control

All support materials included

Sample & pour



Read results immediately!



Product Information

Information Type	Product Details
Method	Flow-through
Time to results	As little as 60 seconds
Storage conditions	35.6–86 °F (2–30 °C)
Shelf life	6 months ¹
Sample type	Fingerstick whole blood, venipuncture whole blood, plasma or serum

Notes

¹Preliminary shelf life. Final shelf life still being established.

Ordering Information

Product

- 90-1062-001 - INSTI® HCV 25 Antibody Test Kits
- 90-1062-002 - INSTI® HCV 50 Antibody Test Kits
- 90-1062-003 - INSTI® HCV 300 Antibody Test Kits
- 90-1105 - INSTI® HCV Bulk Components 50-Pack

Accuracy

Sample Type	Sensitivity	Specificity
Whole blood/Serum/Plasma	100%	99.7%

Manufactured by:



bioLytical Laboratories Inc.

406-13251 Delf Place
Richmond, BC Canada V6V 2A2
1-866-674-6784 | info@biolytical.com

insti.com



65-5001(B) 10/2021